Urine Reagent Strip for

Diagnosis of Veterinary Disease

Available test parameters

Glucose, Protein, pH, Blood, Ketone(Acetoacetic Acid), Specific Gravity, Nitrite, Bilirubin, Urobilinogen, Leukocyte, Microalbumin and Creatinine in Urine



Available diseases diagnosed by the test strip:

-Diabetes
-Kidney(Renal) Problem
-Bacterial Infection
-Urinary Tract Infection
-Liver Disease
-Hepatitis
-Metabolic Disorder
-Etc.



Test Specification	
Product Model Mame	HEALTH MATE™ VET-10, VET-11C, VET-11AC, VET-2AC
Measurement Item	Variable depending on product model name
Required Sample	Approx. 5 mL of Urine
Reaction Time	90~120 sec. for Leukocyte test and 60 sec. for other tests
Reaction Temperature	2~30°C

* HEALTH MATE VET-11AC or VET-2AC strip has Microalbumin and Creatinine tests for diagnosis of early stage of kidney problem with the animals


DFI Co.,Ltd.
388-25, Gomo-ro, Jillye-myeon, Gimhae-si, Gyeongsangnam-do, Korea (South)
Tel: 82-55-3461882 , Fax: 82-55-3461883 , E-mail: sales@cybow.com , Web-site: www.cybow.com

Rev.02/2014-07-01

Veterinary Urinalysis Analyzer

UriDoctor™ VET

Product Features





Specifications

Measurement : Reflectance Photometer
Size(Dimension) : 188 (W) x 74 (D) x 77 (H) mm
Net weight : 0.40 Kg
Power : DC 12V / 3.33A or AAA battery 1.5V x 8 pcs
Throughput : 45(general mode)~120 (quick mode) tests / hour
Interface : USB
Weight : 430g
Memory : 2000 results
Display : 4.3 inch TFT LCD
Software : user program PC s/w is available
Result printing : via PC s/w
Optional item : External printer, Bluetooth

Rev.02/2014-07-01

DFI Co.,Ltd.

388-25, Gomo-ro, Jillye-myeon, Gimhae-si, Gyeongsangnam-do, Korea (South)
Tel: 82-55-3461882 , Fax: 82-55-3461883 , E-mail: sales@cybow.com , Web-site: www.cybow.com